

# ► CUSHIONS



## Water cushion

Ref. W0465 001 001



- By moving under the weight of the body, the water ensures good distribution of pressures.
- Supplied with filling kit.

**Norton score 15 to 20** Pressure ulcer risk **nil to low**



## Foam and fluid gel cushion

Ref. W0470 003 001

- It is composed of 4.5 cm thick 38 kg/m<sup>3</sup> high resiliency foam and a 1 cm thickness of fluid gel. The upper layer of the gel accompanies body movements.
- The foam under the fluid gel ensures an excellent distribution of pressures.
- The cover is impermeable and breathable. It is permeable to water vapor and thus prevents perspiration buildup.
- Cover washable at 90°C, that can be disinfected with any standard detergent.

**WARRANTY**  
**2**  
**YEARS**

**Norton score 13 to 14**

Pressure ulcer risk **low to moderate**

**51**

Ref.	W0465 001 001	W0470 003 001
Cushion weight	(filled) about 9 kg	1.9 kg
Dimensions	43 x 48 cm	40 x 40 x 5.5 cm
Maximum user weight	100 kg	95 kg